

**Chromo  
Pfu DNA Polymerase**



Lot :  
 Expiry Date :  
 Concentration : 1u/μl  
 Supplied with : 2ml of 10X ViBuffer A  
 1ml of 10X ViBuffer S  
 1ml of 50mM MgCl<sub>2</sub>  
 Store at - 20°C

**Product No : PL5205**  
**Quantity : 100u**

info@vivanttechnologies.com



**Description:**

Chromo *Pfu* DNA Polymerase is an extremely thermostable proofreading DNA polymerase. It is suitable for applications requiring high temperature synthesis of DNA. *Pfu* DNA Polymerase catalyzes the polymerization of nucleotides into duplex DNA in the 5' to 3' direction with the presence of Mg<sup>2+</sup>. It exhibits the 3' to 5' proofreading activity, resulting in over 10-fold higher fidelity than possible with *Taq* DNA Polymerases. The enzyme is supplemented with inert color tracer dyes for ease of visualization of the addition of polymerase to the reaction.

**Features:**

- Ultra pure recombinant protein.
- Recommended for use in high-fidelity amplification, amplification of GC-rich sequences or problematic secondary structures, primer extension reactions at elevated temperatures and cloning of blunt-ended amplification products.

**Unit Definition :**

1u is defined as amount of enzyme that required to catalyze the incorporation of 10nmoles of dNTP into acid-insoluble material in 30 minutes at 74°C.

**Reaction Buffer:**

**10X ViBuffer A (without MgCl<sub>2</sub>):**

500mM KCl, 100mM Tris-HCl (pH9.1 at 20°C), 0.1% Triton™X-100. The buffer is optimized for use with 0.1-0.2mM of each dNTP.

**10X ViBuffer S:**


160mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 500mM Tris-HCl (pH 9.2 at 22°C), 17.5mM MgCl<sub>2</sub> and 0.1% Triton™X-100. The buffer is optimized for use with 0.35mM of each dNTP.

**Storage Buffer:**

20mM Tris-HCl (pH 8.0 at 22°C), 100mM KCl, 0.5% Tween™20, 0.5% Nonidet P-40, 0.1mM EDTA, 1mM DTT and 50% glycerol.

**Quality Control:**

All preparations are assayed for contaminating endonuclease, exonuclease, and non-specific DNase activities. Functionally tested in DNA amplification.



**Amplification Using Vivantis *Pfu* DNA Polymerase**  
**Lane M1** : VC 1kb DNA Ladder  
**Lane 0.5 and 1.5kb** : 0.5kb PCR amplification product generated using 0.2mM dNTPs and 2.0u Vivantis *Pfu* DNA Polymerase.  
**Lane 5kb and 8kb** : 5kb and 8kb amplification products generated using 0.25mM dNTPs, 2.5u Vivantis *Pfu* DNA Polymerase and 3% of formamide.  
**Lane M2** : VC Lambda/*Hind*III Marker

0.7% TAE agarose gel

**SUGGESTED INITIAL PCR CONDITIONS FOR VARIOUS PCR PRODUCT SIZES WITH VIVANTIS DNA POLYMERASES (#PL5201 - 06 / #PL5211 - 16)**  
**REACTION MIX (FINAL CONCENTRATION) :**

Primers : 0.5μM Template: Plasmid (0.02-0.2ng) Genomic (0.05-5μg)	Product Size	100bp - 5kb	5kb - 8kb
	dNTP Mix	100μM	200μM
	ViBuffer (1X)	A	S
	Ultrapure DMSO	-	3%
DNA Polymerase	Refer to the below Table (A)		

Top up with sterile dH<sub>2</sub>O to 50 ml.

Product Size	100bp - 5kb	5kb - 8kb
Denaturation	95°C, 2 min	95°C, 2 min
Denaturation	94°C, 20 s	94°C, 12 s
Annealing*	50 - 68°C, 30 s	50 - 68°C, 30 s
Extension / 1kb	72°C, 2 min	72°C, 2 min
Cycles	25 - 35	25 - 35
Final Extension	72°C, 7 min	72°C, 7 min

\* Primer dependent

**TABLE (A) : RECOMMENDED UNITS FOR SPECIFIC VIVANTIS DNA POLYMERASES PER 50μL REACTION VOLUME :**

Product Size	Taq (#PL1201 - 06) <sup>1</sup>	Max Taq (#PL2201 - 06) <sup>2</sup>	Pfu (#PL5201 - 06)	Max Pfu (#PL5211 - 16)
0.1 - 5.0kb	2.0	2.0	2.0	2.0
5.0 - 8.0kb	2.5	2.0	2.5	2.0
8.0 - 20.0kb	2.5	2.0	2.5	2.0
>20.0kb	-	2.0	-	2.0

\* This protocol is subjected to changes depending on the template DNA

<sup>1</sup> Also applied for At Taq (#PL3201-06)  
<sup>2</sup> Also applied for Max At Taq (#PL4201-06)