v i v a n t i s

Nucleotides

Product Datasheet

dNTP Set

Product No: NP2406 Quantity :4 x 0.25ml

Lot : Expiry date :

Concentration: 100mM

Store at -20°C

info@vivantechnologies.com

Description:

dNTP Set consists of 100mM aqueous solution of dATP, dCTP, dGTP and dTTP, each in a separate vial.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V i V a 11 t i S | www.vivantechnologies.com

v i v a n t i S

Nucleotides

Product Datasheet

dNTP Set

Product No: NP2407 Quantity :4 x 1ml

Lot

Expiry date :

Concentration: 100mM

Store at -20°C

info@vivantechnologies.com

Description:

dNTP Set consists of 100mM aqueous solution of dATP, dCTP, dGTP and dTTP, each in a separate vial.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V i V a 11 t i S | www.vivantechnologies.com

v i v a n t i S

Nucleotides

Product Datasheet

dNTP Set

Product No: NP2408 Quantity :4 x 5ml

Lot

Expiry date :

Concentration: 100mM

Store at -20°C

Contact: info@vivantis.com

Description:

dNTP Set consists of 100mM aqueous solution of dATP, dCTP, dGTP and dTTP, each in a separate vial.

Quality Control:

Functionally tested in PCR with *Taq* and *Pfu* DNA Polymerases. Purity of each dNTP >98% by HPLC.



Product Use Limitation

This product is for research purposes and in vitro use only.

V i V a n t i S | www.vivantis.com